

## Isolation of Mouse Primary Sternal Chondrocytes

Karen M. Lyons Lab - UCLA

Modified on 08/12/2011

### Procedure:

- 1) Cut out rib cage from P0-P5 pups and place into complete media (CM: DMEM + 10 % FBS + 1 % Pen/Strep (P/S))
  - a. Euthanize pup by isoflourane inhalation
  - b. Cut out cartilage portion of the rib cage
  - c. Transfer rib cage to a well of a 6-well dish containing PBS +1% P/S
  - d. Repeat a-c for remaining pups
  - e. Replace PBS w/ CM and incubate rib cage overnight in 37 °C/5% CO<sub>2</sub> or can continue to step 2
- 2) Digest rib cage in 12-well dish or 15 ml conical tube:
  - a. Make protease (Sigma P6911) digest mix:
    - i. Per rib cage (2 ml):
      1. 2 mg/ml protease in PBS + 1% P/S
      - 2. Filter sterilize w/ 0.45 um filter**
  - b. Add 2 ml protease digest mix to each rib cage
    - i. Incubate at 37 °C incubator for 15 minutes
  - c. Rinse ribs 1X w/ PBS + 1% P/S
  - d. Make 3 mg/ml collagenase (Sigma C6885) digest mix:
    - i. Per rib cage (2 ml):
      1. 3 mg/ml collagenase in DMEM + 1% P/S
      - 2. Filter sterilize w/ 0.45 um filter**
  - e. Add 2 ml of 3 mg/ml collagenase digest mix to each rib cage
    - i. Incubate at 37 °C incubator for 15 minutes
  - f. Rinse each rib 3X w/ PBS + 1% P/S in 15 ml conical tube
    - i. Let cartilage settle to the bottom of the tube
    - ii. Aspirate PBS and soft tissue on top layer
    - iii. Transfer cartilage to a well of a new 12-well plate
  - g. Make 0.3 mg/ml collagenase digest mix:
    - i. Per rib cage (3 ml for 12-well plate):
      1. 0.3 mg/ml collagenase in DMEM + 1% P/S
      - 2. Filter sterilize w/ 0.45 um filter**
  - h. Add 3 ml of 0.3 mg/ml collagenase digest mix to each rib cage
    - i. Incubate overnight at 37 °C/5% CO<sub>2</sub>
  - i. Add 3 ml of CM to each rib cage
  - j. Filter through 70 µm cell strainer onto a clean 50 ml conical tube
    - i. Centrifuge at 1000 rpm for 5 min
  - k. Wash 2X w/ CM
    - i. Centrifuge at 1000 rpm for 5 min between washes
  - l. Seed cells onto a well of a 6-well plate at  $1 \times 10^6$  cells/well or a 60 mm dish at  $5 \times 10^6$  cells/well in chondrogenic media (CM + 50 µg/ml ascorbic acid + 10 mM glycerol-2-phosphate)
  - m. Change media the following day
  - n. Expand to 100% confluency