

# Strategy for Determining the Activity of a Gene's Promoter

Performed by HC70AL students

## Arabidopsis Genomic DNA

- ↓ **Amplify** the promoter by PCR using High Fidelity DNA Polymerase and a modified FW primer beginning with CACC

## Amplified product beginning with CACC

Linearized pENTR vector bound by Topoisomerase, also containing GTGG overhang

- ↓ CACC in PCR product **anneals** to complementary GTGG in pENTR vector
- ↓ Topoisomerase **seals** the phosphodiester bond and is released from the plasmid

## Population of recombinant and non-recombinant plasmid

- ↓ **Transform** *E. coli* cells using CaCl
- ↓ **Screen** for colonies containing recombinant plasmid by isolating plasmid DNA and performing a restriction digest and gel electrophoresis
- ↓ **Sequence** the promoter to make sure there are no mutations

## Promoter cloned into pENTR vector

pBGWFS7 T-DNA vector (T-DNA contains toxic gene upstream of *GFP* and *GUS* reporter genes, BastaR gene, and left and right border sequences)

- ↓ Homologous recombination between the *att* sequences flanking the promoter in pENTR and flanking the toxic gene in pBGWFS7

## T-DNA plasmid (T-DNA contains promoter upstream of *GFP* and *GUS* reporter genes, BastaR gene, and left and right border sequences)

- ↓ **Transform** *E. coli* cells using CaCl
- ↓ **Screen** for colonies containing recombinant plasmid by isolating plasmid DNA and performing a restriction digest and gel electrophoresis

## Promoter cloned into T-DNA plasmid

- ↓ **Transform** *Agrobacterium tumefaciens* cells containing the Ti Helper plasmid, which encodes genes with T-DNA transfer functions
- ↓ **Transform** *Arabidopsis* plant with *Agrobacterium* cells containing Ti Helper plasmid and the T-DNA plasmid
- ↓ **Select** for transformed *Arabidopsis* plants by applying the herbicide Basta

Transformed *Arabidopsis* plants containing the promoter upstream of *GFP* and *GUS* genes